525 Series  
Serie de 525  
Série 525

## **EN** DeVilbiss® 5 Liter Oxygen Concentrator Instruction Guide

**WARNING**– READ INSTRUCTION GUIDE BEFORE OPERATING THIS EQUIPMENT.

**CAUTION**– *Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.*

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**CAUTION**– *Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.*

**WARNING**

**Under certain circumstances, oxygen therapy can be hazardous. Seek medical advice before using an oxygen concentrator.**

**Physician Information**

Physician Name: \_\_\_\_\_

Telephone: \_\_\_\_\_

Address: \_\_\_\_\_

**Prescription Information**

Name: \_\_\_\_\_

Oxygen liters per minute

    at rest: \_\_\_\_\_ during activity: \_\_\_\_\_ other: \_\_\_\_\_

Oxygen use per day

    Hours: \_\_\_\_\_ Minutes: \_\_\_\_\_

    Comments: \_\_\_\_\_


**DeVilbiss 5-Liter Oxygen Concentrator w/OSD Serial Number:** \_\_\_\_\_

**DeVilbiss Equipment Provider Information**

Set-Up Person: \_\_\_\_\_






This instruction guide was reviewed with me, and I have been instructed on the safe use and care of the DeVilbiss Oxygen Concentrator.

Signature: \_\_\_\_\_ Date: \_\_\_\_\_

**DeVilbiss 5-Liter Series**

## IEC SYMBOLS

	It is mandatory to read and understand the operating instructions prior to use.		Double Insulated		Danger - No smoking near patient or device.
	Keep unit dry.		Catalog Number		Do NOT use Oil, Grease or Lubricants
	Electric Shock Hazard. Cabinet to be removed by authorized personnel only.		Serial Number		Do not use near heat or open flames
	Off On		External Power		Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
	Reset		Normal Oxygen		Manufacturer
	Alternating Current		Low Oxygen		Maximum recommended flow rate: 5LPM
	Type B applied part		Service Required		
	This device contains electrical and/or electronic equipment that must be recycled per EU Directive 2012/19/EU- Waste Electrical and Electronic Equipment (WEEE)				

## IMPORTANT SAFEGUARDS

Read this entire guide before using your DeVilbiss concentrator. Important safeguards are indicated throughout this guide. Pay special attention to all safety information.

### READ ALL INSTRUCTIONS BEFORE USING.

#### WARNING:

*A warning indicates the possibility of injury to the user or the operator.*

- Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy.
- Keep the oxygen concentrator and cannula at least 2 m (6.5 feet) from hot, sparking objects or naked sources of flame.
- To ensure you receive the therapeutic amount of oxygen delivery according to your medical condition, 525 Oxygen Concentrator must:
  - be used only after one or more settings have been individually determined or prescribed for you at your specific activity levels.
  - be used with the specific combination of parts and accessories that are in line with the specification of the concentrator manufacturer and that were used while your settings were determined.
- For your safety, the oxygen concentrator must be used according to the prescription determined by your physician.
- Your delivery settings of the oxygen concentrator should be periodically reassessed for the effectiveness of therapy.
- Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- This device contains electrical and/or electronic equipment. Follow local governing ordinances and recycling plans regarding disposal of device components.
- Installation of 515LF-607 low output flow meter package will cause the low flow alarm to not work.
- Use only water-based lotions or salves that are oxygen-compatible before and during oxygen therapy. Never use petroleum or oil-based lotions or salves to avoid the risk of fire and burns.
- Do not lubricate fittings, connections, tubing or other accessories of the oxygen concentrator to avoid the risk of fire and burns.
- Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns.
- Use of this device at an altitude above 4000 meters or above a temperature of 95°F (35°C) or greater than 95% relative humidity may affect the flow rate and the percentage of oxygen and consequently the quality of the therapy. Refer to specifications for details regarding parameters tested.
- Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula or mask on bed coverings or chair cushions if the oxygen concentrator is turned on but not in use. The oxygen will make the materials flammable. Turn the oxygen concentrator off when not in use to prevent oxygen enrichment.
- If you feel discomfort or are experiencing a medical emergency while undergoing oxygen therapy, seek medical assistance immediately to avoid harm.
- Geriatric, pediatric or any other patient unable to communicate discomfort can require additional monitoring and/or a distributed alarm system to convey the information about the discomfort and/or the medical urgency to the responsible caregiver to avoid harm.
- Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking within the same room where the oxygen concentrator or any oxygen carrying accessories are located. If you intend to smoke, you must always turn the oxygen concentrator off, remove the cannula and leave the room where either the cannula or mask or the oxygen concentrator is located. If unable to leave the room, you must wait 10 minutes after you have turned off the oxygen concentrator before smoking.
- Open flames during oxygen therapy are dangerous and are likely to result in fire or death. Do not allow open flames within 2 m (6.5 feet) of the oxygen concentrator or any oxygen carrying accessories.
- Locate oxygen tubing and power supply cords to prevent tripping hazards and reduce the possibility of entanglement or strangulation.
- The Device is classified as IPX0 which means it does not provide protection from water ingress. Keep the device dry.
- Do not service or clean this device while in use with a Patient.
- Under certain circumstances, oxygen therapy can be hazardous. Seek medical advice before using an oxygen concentrator.
- It is recommended that the homecare provider lock the flow control knob to prevent inadvertent adjustment. A flow setting other than prescribed may affect the patient therapy.
- When device is used under extreme operating conditions, the temperature near the exhaust vents on the bottom of the unit may reach 145°F (63°C). Keep body parts a minimum of 32" away from this area.
- Improper use of the power cord and plugs can cause a burn, fire or other electric shock hazards. Do not use the unit if the power cord is damaged.
- The DeVilbiss oxygen concentrators manufactured since July 1, 2012 are equipped with a fire mitigating outlet fitting that prevents propagation of fire into the unit. DeVilbiss oxygen concentrators manufactured prior to July 1, 2012 must be retrofitted by your provider with the fire mitigating outlet fitting.
- In order to prevent a fire propagating from the patient through the cannula towards the unit, a means of protection should be located as close to the patient as practicable. Please contact your provider for this means of protection.
- Before attempting any cleaning procedures, turn the unit "Off."

## NO SMOKING

- To avoid electric shock, do not remove the concentrator cabinet. The cabinet should only be removed by a qualified DeVilbiss technician. Do not apply liquid directly to the cabinet or utilize any petroleum-based solvents or cleaning agents.
- Use of harsh chemicals (including alcohol) is not recommended. If bactericidal cleaning is required, a non-alcohol based product should be used to avoid inadvertent damage.
- Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the Electromagnetic Compatibility [EMC] information provided in the accompanying documents.
- Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.
- The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

**CAUTION:**

*A Caution indicates the possibility of damage to the device.*

- It is very important to follow your oxygen prescription. Do not increase or decrease the flow of oxygen – consult your physician.
- To prevent product damage, do not attempt to operate the unit without the air filter or while the filter is still damp.

# SAVE THESE INSTRUCTIONS.

## INTRODUCTION

This instruction guide will acquaint you with the DeVilbiss oxygen concentrator. Make sure that you read and understand this guide before operating your unit. Important safeguards are indicated throughout this guide. Pay special attention to all safety information. Contact your DeVilbiss equipment provider should you have any questions.

### Intended Use

The DeVilbiss 5 Liter Oxygen Concentrator intended use is to provide supplemental low flow oxygen therapy for patients suffering from COPD, cardiovascular disease and lung disorders. The oxygen concentrator is used in home type environments, homes, nursing homes, patient care facilities, etc.

### Indications For Use

The DeVilbiss Oxygen Concentrator is intended for use as an oxygen concentrator to provide supplemental low flow oxygen therapy in the home, nursing homes, patient care facilities, etc.

**NOTE**– The device is not intended to be life supporting or life sustaining.

### Why Your Physician Prescribed Supplemental Oxygen

Today, many people suffer from heart, lung and other respiratory diseases. Many of these people can benefit from supplemental oxygen therapy. Your body requires a steady supply of oxygen to function properly. Your physician prescribed supplemental oxygen for you because you are not getting enough oxygen from room air alone. Supplemental oxygen will increase the amount of oxygen that your body receives.

Supplemental oxygen is not addictive. Your physician prescribed a specific oxygen flow to improve symptoms such as headaches, drowsiness, confusion, fatigue or increased irritability. If these symptoms persist after you begin your supplemental oxygen program, consult your physician.

The oxygen delivery setting has to be determined for each patient individually with the configuration of the equipment to be used, including accessories.

The proper placement and positioning of the prongs of the nasal cannula in the nose is critical to the amount of oxygen delivered to the respiratory system of the patient.

Your Delivery settings of the oxygen concentrator should be periodically reassessed for the effectiveness of therapy.

### How Your DeVilbiss Oxygen Concentrator Works


Oxygen concentrators are the most reliable, efficient and convenient source of supplemental oxygen available today. The oxygen concentrator is electrically operated. The unit separates oxygen from room air which allows high-purity supplemental oxygen to be delivered to you through the oxygen outlet. Although the concentrator filters the oxygen in a room, it will not affect the normal amount of oxygen in your room.

## IMPORTANT PARTS OF YOUR CONCENTRATOR

Please take time to familiarize yourself with your DeVilbiss oxygen concentrator before operating.


### Front View (Figure A)

1. Operating instructions
2. Green Power light – illuminates when your concentrator is operating
3. Power Switch  
| = ON  
O = OFF
4. Flow meter knob
5. Flow meter
6. Circuit breaker – resets the unit after electrical overload shutdown
7. Oxygen outlet – oxygen is dispersed through this port
8. Normal Oxygen (green) light (see page 6)
9. Low Oxygen (yellow) light (see page 6)
10. Red Service Required light – when illuminated contact your DeVilbiss provider

### Back View (Figure B)

11. Handgrip
12. Exhaust
13. Power cord and/or IEC power connector.
14. Line cord strap
15. Air filter – prevents dirt, dust and lint from entering your unit.
16. Auxiliary Oxygen Port (N/A for 525PS): Your concentrator is equipped with an auxiliary oxygen port that can be used to fill oxygen cylinders with an FDA-cleared cylinder filling device that is designed to use oxygen from a concentrator to fill a cylinder. The port is only for use with FDA-cleared filling devices with compatible oxygen input specifications. Refer to the cylinder filling device instruction guide for the oxygen input/output specifications, connection and operating instructions.

## Accessories

Transfiller Caddy ..... DeVilbiss 525DD-650  
Oxygen Outlet Connector - Plastic - 1/pack..... CN100  
Bubble Humidifier ..... Salter Labs 7600 or equivalent

There are many types of humidifiers, oxygen tubing and cannulas/masks that can be used with this device. Certain humidifiers and accessories may impair the device's performance. A mask or any nasal cannula can be used with continuous flow delivery and may be sized according to your prescription as recommended by your homecare provider who should also give you advice on the proper usage, maintenance and cleaning.

### WARNING:

When using the Transfiller Caddy with a Transfill device, always keep the system on a flat surface. Disassemble the system prior to moving.

**NOTE**– The bubble humidifier should be supplied with a permanent fire stop device. If a bubble humidifier needs to be used without a permanent fire stop device, a secondary fire stop device must be used and placed as close to the humidifier as possible. Failing to do so could increase the risk of fire.

**NOTE**– A maximum of 50 feet (15 meters) of crush-proof oxygen tubing, plus 7 feet (2.1 meters) of cannula, plus a bubble humidifier is allowed between the concentrator and the patient.

**NOTE**– The oxygen supply accessory (patient tubing) shall be equipped with a means that, in case of fire, stops the delivery of oxygen to the patient. This means of protection should be located as close to the patient as practicable.

**NOTE**– Your healthcare provider should verify the compatibility of the oxygen concentrator and all of the parts used to connect to the patient before use.

## SETTING UP YOUR OXYGEN CONCENTRATOR

1. Position your unit near an electrical outlet in the room where you spend most of your time.

**NOTE**– Do not connect to an electrical outlet controlled by a wall switch.

### WARNING

Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy. Keep the oxygen concentrator and cannula at least 2 m (6.5 feet) from hot, sparking objects or naked sources of flame.

2. Position your unit on a flat surface at least 6 inches (16 cm) from walls, draperies or any other objects that might prevent the proper flow of air in and out of your oxygen concentrator. The oxygen concentrator should be located so as to avoid pollutants or fumes.  
**NOTE**– To move the unit, firmly grasp the handle located on the top of the unit, rolling and/or lifting the unit over pathway obstacles.
3. Before operating your unit, always check to be sure the air filter (located on the back of your unit) is clean. Proper cleaning of this filter is discussed in the Caring For Your Concentrator section on page 6.
4. Attach the appropriate oxygen accessories to the oxygen outlet.

#### Oxygen Tubing Connection:

- a. Thread the oxygen outlet connector onto the oxygen outlet.
- b. Attach the oxygen tubing directly to the connector (Figure 1).

#### Oxygen Tubing Connection With Humidification:

If your physician has prescribed an oxygen humidifier as part of your therapy, follow these steps (If using a prefill, go to step b.):

- a. Fill the humidifier bottle with distilled water. Do not overfill.
- b. Thread the wing nut located on the top of the humidifier bottle to the oxygen outlet so that it is suspended (Figure 2). Make sure it is securely tightened.
- c. Attach the oxygen tubing directly to the humidifier bottle outlet fitting (Figure 3).

**NOTE**– Your physician has prescribed either a nasal cannula or face mask. In most cases, they are already attached to the oxygen tubing. If not, follow the manufacturer's instructions for attachment.

5. Remove the power cord completely from the line cord strap. Make sure the power switch is in the "Off" position and insert the plug into the wall outlet. The unit is double insulated to guard against electric shock.

**NOTE**– (only 115 volt units) The plug on the DeVilbiss oxygen concentrator has one blade wider than the other. To reduce the risk of electric shock, this plug is intended to fit in a wall outlet only one way. Do not attempt to defeat this safety feature.

**NOTE**– To check your oxygen concentrator and accessories for proper operation; 1. Check the output flow by placing the end of the nasal cannula under the surface of a half-full cup of water and look for bubbles. 2. Check the system for leaks by bending the nasal prongs over and squeeze tight to stop the flow of oxygen. Look at the flow meter to see that the indicator ball on the flow meter drops to zero. If the indicator ball does not drop to zero, check all connections for possible leaks. Parts to check for leaks are: tubing connections, humidifier bottle and other accessories like firebreaks. Repeat these steps until the flow meter ball drops to zero. Contact your provider or service supplier immediately if you encounter any problems.

### WARNING

Improper use of the power cord and plugs can cause a burn, fire or other electric shock hazards. Do not use the unit if the power cord is damaged.

## OPERATING YOUR DEVILBISS OXYGEN CONCENTRATOR

### WARNING

Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy. Keep the oxygen concentrator and cannula at least 2 m (6.5 feet) from hot, sparking objects or naked sources of flame.

The DeVilbiss oxygen concentrators manufactured since July 1, 2012 are equipped with a fire mitigating outlet fitting that prevents propagation of fire into the unit. DeVilbiss oxygen concentrators manufactured prior to July 1, 2012 must be retrofitted by your provider with the fire mitigating outlet fitting.

In order to prevent a fire propagating from the patient through the cannula towards the unit, a means of protection should be located as close to the patient as

**NO SMOKING**

practicable. Please contact your provider for this means of protection.

1. Press the power switch to the "On" position. When the unit is turned "On," all four lights (Power, Service Required, Low Oxygen and Normal Oxygen) on the front panel will illuminate briefly, and an audible signal will briefly alarm. After a few seconds, only the Power and Normal Oxygen lights will remain lit.

**NOTE**– DeVilbiss recommends for optimal service life that the DeVilbiss Oxygen Concentrator be operated for at least 30 minutes after it is powered on. Shorter periods of operation, operating in extreme temperature/humidity conditions or in the presence of contaminants and/or handling and storage conditions outside those specified, may affect the long term reliable operation of the product.

### **WARNING**

**Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula or mask on bed coverings or chair cushions if the oxygen concentrator is turned on but not in use. The oxygen will make the materials flammable. Turn the oxygen concentrator off when not in use to prevent oxygen enrichment.**

**NOTE**– If the Service Required light illuminates, and the audible signal alarms but the unit is not operating, there is no power to the unit. Refer to the Minor Troubleshooting chart on page 7 and contact your DeVilbiss provider if necessary.

**NOTE**– If an audible low-frequency vibration sound is detected, the unit is not operating properly. Refer to the Minor Troubleshooting chart on page 7 and contact your DeVilbiss provider if necessary.

2. Check the flow meter to make sure that the flow meter ball is centered on the line next to the prescribed number of your flow rate.

**CAUTION**– It is very important to follow your oxygen prescription. Do not increase or decrease the flow of oxygen – consult your physician.

**NOTE**– Your DeVilbiss provider may have preset the flow meter so that it can not be adjusted.

**NOTE**– If the flow meter knob is turned clockwise, the flow decreases (and eventually will shut off the oxygen flow). If the knob is turned counterclockwise, the flow increases.

**NOTE**– For prescriptions of 5 LPM, be sure the ball is centered on the 5 liter line. The ball should not touch the red line. Setting the flow higher than 5 may cause the oxygen purity level to drop.

**NOTE**– The low-flow alarm may activate if the flow meter ball is set at or below 0.5 lpm. The unit will continue to run; however, the Service Required light will come on accompanied by an audible alarm. Adjust the flow meter to your prescribed flow.

**NOTE**– Installation of 515LF-607 low output flow meter package will cause the low flow alarm to not work.

3. Your DeVilbiss concentrator is now ready for use. Properly position the cannula or mask (Figure 4) and allow 20 minutes for the oxygen concentrator to reach stated performance.

### **DeVilbiss OSD® Operation**

The OSD (Oxygen Sensing Device) is a device within your concentrator that monitors the oxygen produced by your unit. When the unit is turned "On," all four lights (Power, Service Required, Low Oxygen and Normal Oxygen) on the front panel will illuminate briefly. After a few seconds, only the Power and Normal Oxygen lights will remain lit.

The OSD lights on the front panel are defined as follows:

- Green Normal Oxygen light–acceptable oxygen level.
- Yellow Low Oxygen light– below an acceptable oxygen level.

If the oxygen purity falls below the acceptable level, the green Normal Oxygen light will shut off, and the yellow Low Oxygen light will illuminate. Switch to your reserve oxygen system. Refer to the Minor Troubleshooting section in this guide on page 7 and contact your DeVilbiss provider.

As an added safety feature should the oxygen purity continue to drop, an intermittent audible signal will sound. Contact your DeVilbiss provider immediately. Do not attempt any other maintenance.

## **RESERVE OXYGEN SYSTEM**

As a precaution, your DeVilbiss provider may supply you with a reserve oxygen system. If your unit loses electrical power or fails to operate correctly, the Patient Alert System will sound to signal you to switch to your reserve oxygen system (if provided) and contact your DeVilbiss provider.

## **CARING FOR YOUR DEVILBISS OXYGEN CONCENTRATOR**

DeVilbiss recommends using only original DeVilbiss parts and filters in order to guarantee reliable operation of the product.

### **WARNING**

**Do NOT use lubricants, oils or grease.**

**Before attempting any cleaning procedures, turn the unit "Off."**

### **Cannula/Mask, Tubing and Humidifier Bottle**

Clean and replace the cannula/mask, tubing and humidifier bottle according to the manufacturer's instructions.

### **Air Filter and Oxygen Outlet Connector**

The air filter and connector should be cleaned at least once a week. To clean, follow these steps:

1. Remove the air filter located on the back of the unit. Remove the oxygen outlet connector (if used).
2. Wash in a solution of warm water and dishwashing detergent (Figure 5).
3. Rinse thoroughly with warm tap water and towel dry. The filter should be completely dry before reinstalling.

**CAUTION**– To prevent product damage, do not attempt to operate the unit without the air filter or while the filter is still damp.

### **Exterior Cabinet**

Clean the concentrator exterior cabinet weekly by using a damp cloth and wiping dry.

### **WARNING**

**To avoid electric shock, do not remove the concentrator cabinet. The cabinet should only be removed by a qualified DeVilbiss technician. Do not apply liquid directly to the cabinet or utilize any petroleum-based solvents or cleaning agents.**

**Use of harsh chemicals (including alcohol) is not recommended. If bactericidal cleaning is required, a non-alcohol based product should be used to avoid inadvertent damage.**

## WARRANTY

DeVilbiss provides warranty for this device. Use only DeVilbiss approved original spare parts for maintenance and repair to retain warranty. Using unapproved or non-original spare parts would void warranty.

## PROVIDER'S NOTES - Cleaning and Disinfection When There is a Patient Change

When medical devices have already been used with a patient, contamination with human pathogenic germs should be assumed (unless there is evidence to the contrary). The next patient, user or third party should be protected by appropriate handling and preparation.

Therefore, when there is a patient change, people must be protected during the transport and handling of the device, and the device must be fully processed, i.e., cleaned and disinfected, by suitably trained personnel before reuse to protect the next patient. The complete processing may only be done by the manufacturer or by a qualified DeVilbiss provider/service technician.

**NOTE**– If the following described complete processing of the concentrator by a qualified DeVilbiss provider/technician is not possible, the device must not be used by another patient!

DeVilbiss Healthcare recommends that at least the following procedures be carried out by the manufacturer or a qualified third party between uses by different patients.

**NOTE**–If preventive maintenance is due at this time, these procedures should be carried out in addition to the servicing procedures.

1. Dispose of all accessory components that are not suitable for reuse, i.e., particularly the oxygen tubing, the nasal cannula/mask, oxygen outlet connector and humidifier bottle.
2. **The concentrator must be disconnected from the power supply for this step:** Open the concentrator and remove all dust deposits inside the cabinet with an appropriate vacuum cleaner.
3. Clean and disinfect the power cord and all internal and external parts with a suitable disinfecting agent, e.g., Microbac Forte or Terralin®.
4. Check the cord, the plug on the back of the device, the power switch, the fuse holder and the indicator light for possible damage.
5. Replace all damaged or worn components.
6. Replace the cabinet air filter on the back of the device.
7. Check the oxygen concentration. If the device is within specifications, the extended life intake bacteria filter does not need to be replaced between patients. If the concentration is not within specifications, the provider should refer to the service manual section on Troubleshooting.

**NOTE**– There is no portion of the gas pathways through the concentrator that could be contaminated with body fluids under normal conditions.

The device patient connection may unintentionally become contaminated with expired gases for a single fault condition i.e., a hose internal to the device becomes disconnected. This condition will cause no flow out of the device and/or an alarm condition. Should this occur refer to the Service Manual for additional instructions.

## TROUBLESHOOTING

The following troubleshooting chart will help you analyze and correct minor oxygen concentrator malfunctions. If the suggested procedures do not help, switch to your reserve oxygen system and call your DeVilbiss homecare provider. Do not attempt any other maintenance.

### WARNING

To avoid electric shock hazard, do not remove the cabinet. The cabinet should only be removed by a qualified DeVilbiss homecare technician.

### Minor Troubleshooting Chart

SYMPTOM	POSSIBLE CAUSE	REMEDY
A. Unit does not operate. Power light is off when the power switch is "On." Audible alert is pulsing and Service Required light is flashing.	1. Power cord not properly inserted into wall outlet.	1. Check power cord connection at the wall outlet. On 230 volt units, also check the mains connection on the back of the unit.
	2. No power at wall outlet.	2. Check your home circuit breaker and reset if necessary. Use a different wall outlet if the situation occurs again.
	3. Oxygen concentrator circuit breaker activated.	3. Press the concentrator circuit breaker reset button located below the power switch. Use a different wall outlet if the situation occurs again. If the above remedies do not work, contact your DeVilbiss provider.
B. Unit operates, the Power light is on when the power switch is "On." Red Service Required light is illuminated. Audible alert may be sounding.	1. Air filter is blocked.	1. Check the air filter. If the filter is dirty, wash it following the cleaning instructions on page 6.
	2. Exhaust is blocked.	2. Check the exhaust area and make sure there is nothing restricting the unit exhaust.
	3. Blocked or defective cannula, face mask or oxygen tubing.	3. Detach cannula or face mask. If proper flow is restored, clean or replace if necessary. Disconnect the oxygen tubing at the oxygen outlet. If proper flow is restored, check oxygen tubing for obstructions or kinks. Replace if necessary.
	4. Blocked or defective humidifier bottle.	4. Detach the humidifier from the oxygen outlet. If proper flow is obtained, clean or replace humidifier.
	5. Flow meter set too low.	5. Set flow meter to prescribed flow rate. If the above remedies do not work, contact your DeVilbiss provider.
C. Unit operates, the power light is on when power switch is "on," audible low-frequency vibration sound is detected.		1. Turn your unit "Off." Switch to your reserve oxygen system and contact your DeVilbiss provider immediately.
D. Both the green Normal Oxygen and the yellow Low Oxygen lights are either on or off.	1. OSD malfunction.	1. Contact your DeVilbiss provider.
E. Yellow Low Oxygen light is on or the yellow Low Oxygen light is on, and the intermittent audible signal is sounding.	1. Flow meter is not properly set.	1. Ensure the flow meter is properly set to the prescribed number. (The maximum flow meter setting is 3 LPM when an oxygen bottle is being filled with oxygen from the auxiliary port.)
	2. Air filter is blocked.	2. Check the air filter. If the filter is dirty, wash it following the cleaning instructions on page 6.
	3. Exhaust is blocked.	3. Check the exhaust area and make sure there is nothing restricting the unit exhaust. If the above remedies do not work, contact your DeVilbiss provider.

SYMPTOM	POSSIBLE CAUSE	REMEDY
F. Red Service Required light is on, and an intermittent audible signal is sounding.	1. Flow meter is not properly set.	1. Ensure the flow meter is properly set to the prescribed number. (The maximum flow meter setting is 3 LPM when an oxygen bottle is being filled with oxygen from the auxiliary port.)
	2. Air filter is blocked.	2. Check the air filter. If the filter is dirty, wash it following the cleaning instructions on page 6.
	3. Exhaust is blocked.	3. Check the exhaust area and make sure there is nothing restricting the unit exhaust. If the above remedies do not work, contact your DeVilbiss provider.
G. If any other problems occur with your oxygen concentrator.		1. Turn your unit "Off." Switch to your reserve oxygen system and contact your DeVilbiss provider immediately.

## RETURN AND DISPOSAL

This device may not be disposed of with household waste. After use of the device, please return the device to the provider for disposal. This device contains electrical and/or electronic components that must be recycled per EU Directive 2012/19/EU-Waste Electrical and Electronic Equipment (WEEE). Non-infectious used accessories (e.g. nasal cannula) can be disposed of as residential waste. The disposal of infectious accessories (e.g. nasal cannula from an infected user) must be made via an approved waste disposal company. Names and addresses can be obtained from the local municipality.

## SPECIFICATIONS

DEVILBISS 5-LITER SERIES	
<b>Catalog Number</b>	<b>525DS, 525DS-Q</b>
<b>Delivery Rate</b> (Lower delivery rates available for low flow applications)	0.5 to 5 LPM
<b>Maximum Recommended Flow</b> (@ nominal outlet pressures of zero and 7 kPa)**	5 LPM
<b>Outlet Pressure</b>	8.5 ± 0.5 psig (58.6 ± 3.5 kPa)
<b>Auxiliary Oxygen Port **</b>	Outlet Pressure: <15 psi Outlet Flow: 2 LPM
<b>Electrical Rating</b>	115 V, 60 Hz, 3.3 Amp
<b>Operating Voltage Range</b>	97-127 V~, 60 Hz
<b>Oxygen Percentage</b>	1-5 LPM=87%-96%
<b>Operating Altitude</b>	
(tested at 70°F {21°C} only) 0-1500 M (0-4921 ft)	Across the voltage range: No degradation of performance
1500-4000 M (4921-13123 ft)	Tested at nominal voltage only: No degradation of performance
<b>Operating Environment Range*</b>	
41°F (5°C) to 95°F (35°C), humidity range of 10% to 95%	No degradation in performance across the operating voltage range.
<b>Power Consumption</b>	310 Watts Average 275 Watts @ 1.2 LPM & below
<b>Weight</b>	36 lbs. (16.3 Kilograms)
<b>Safe Working Load</b>	53 lbs. (24 Kilograms)
<b>Sound Level</b> (ISO 8359:1996)	48 dBA (525DS) 46 dBA (525DS-Q)
<b>Dimensions</b>	24.5"H x 13.5"W x 12"D (62.2 x 34.2 x 30.4 cm)
<b>Pressure Relief Valve</b>	45 psig±5psig (310 kPa±34.5 kPa)
<b>Operating System</b>	Time Cycle / Pressure Swing
<b>The visible "low oxygen" indicator will activate at the following level</b>	84% ± 2% (The audible alarm will alert at approximately 75%. At less than 60%, the red "service required" light will activate.)
<b>Storage Conditions</b>	-40°F (-40°C) to 158°F (70°C), humidity range of 10% to 100%, including condensation
<b>Equipment Class and Type</b>	<input type="checkbox"/> Class II Equipment Double Insulated;  Type B Applied Part Ordinary Equipment, IPX0
<b>Approval Body and Safety Standard</b>	TUV IEC 60601-1 3rd Edition
<b>EMC Compliance To</b>	EN60601-1-2

\*NOTE– The OSD performance at 41°F (5°C) to 104°F (40°C), 95% R.H. through voltage range on the 525DS verified at 670m.

\*\*The maximum recommended flow is 3 LPM when an oxygen bottle is being filled with oxygen from the auxiliary oxygen port.

Specifications subject to change without notice.

## DEVILBISS GUIDANCE AND MANUFACTURER'S DECLARATION

### WARNING

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the Electromagnetic Compatibility [EMC] information provided in the accompanying documents.

Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

The equipment or system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

**NOTE**– The EMC tables and other guidelines provide information to the customer or user that is essential in determining the suitability of the Equipment or System for the Electromagnetic Environment of use, and in managing the Electromagnetic Environment of use to permit the Equipment or System to perform its intended use without disturbing other Equipment and Systems or non-medical electrical equipment.

### Guidance and Manufacturer's Declaration – Electromagnetic Emissions

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1	This device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class B	
Voltage fluctuations / flicker emissions	Complies	

### Guidance and Manufacturer's Declaration – Electromagnetic Immunity

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	Complies	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2.5GHz	Complies	Field strengths outside the shielded location from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than 3 V/m. Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted RF IEC 61000-4-6	3 Vrms 150kHz to 80MHz	Complies	
Electrical fast transient IEC 61000-4-4	±2kV power line ±1kV I/O lines	Complies	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential ±2kV common	Complies	
Power frequency magnetic field IEC 61000-4-8	3 A/m	Complies	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interrupts and voltage variations on power supply input lines IEC 61000-4-11	>95% dip 0.5 cycle 60% dip 5 cycles 70% dip 25 cycles 95% dip 5 secs.	Complies	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or battery.

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DeVilbiss® and OSD® are registered trademarks of DeVilbiss Healthcare.
© 2017 DeVilbiss Healthcare LLC. 11.17 All Rights Reserved.
**DeVilbiss Healthcare Ltd**
**DeVilbiss Healthcare Pty. Limited**
**DeVilbiss Healthcare SAS**
**DeVilbiss Healthcare GmbH**