GUARDIAN[®]

ROLLATOR ROLLING WALKER WITH 8" WHEELS

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GUARDIAN®

Contents

Safety Instructions	3
Assembly Instructions	4-5
Operating Instructions	6-7
Maintenance	7

Safety Instructions

- **DO NOT** use this product without first reading and understanding the instructions contained in this booklet. If you are unable to understand the warnings, cautions, or instructions, contact a healthcare professional, dealer, or technical personnel before use—otherwise serious bodily injury and damage may occur.
- Consult with your doctor, physical, or occupational therapist or other qualified professional to assure the proper height adjustment.
- **DO NOT** use the rollator as a wheelchair
- **DO NOT** attempt to move the rollator while you or anyone is sitting on the seat.
- **DO NOT** self-propel the rollator while seated.
- Serious injury to you and/or damage to the rollator walker frame or wheels may result from improper use.
- **DO NOT** use the seat to transport people or objects. **DO NOT** use the seat to carry or move anything.
- **DO NOT** sit on the rollator while it is being loaded into a vehicle or while it is being transported in a vehicle.
- **DO NOT** exceed the maximum weight capacity (including bag contents) or serious injury could result: **Maximum weight capacity: 300 lbs (136 kg)**
- Before each use make sure the rollator frame is fully opened.
- Check the brakes for proper operation. Make sure that all parts are secure and that moving parts are in good working order. Wheels should roll and not wobble.
- All wheels must be in contact with the floor and/or ground at **ALL TIMES**. This will ensure the rollator is properly balanced.
- The brakes **MUST** be in the locked position **BEFORE** sitting on the seat.
- **DO NOT** sit on the walker if ground is sloped or uneven.
- **DO NOT** use the rollator to walk backwards, down gradients or to climb stairs, curbs, or to go over obstacles. Serious risk of fall or injury may occur.
- **DO NOT** perform any adjustments to the rollator while it is in use.
- Only use accessories and spare parts authorized by Medline Industries, Inc.
- **DO NOT** push rollator while brakes are enabled. Brakes are used to keep the wheels from rolling. But it may not prevent the rollator from moving if it is pushed.
- **DO NOT** use the rollator to provide assistance while rising from a lower position. If you use it to pull yourself up to a standing position, the rollator may move even if the brakes are in the locked position and this could result in a fall.
- Inspect the handle bar grips prior to use. Assure the plastic parts are still functional if they are exposed to extreme temperatures (above 100°F (38°C) or below 32°F (0°C) high humidity and/or become wet.
- **DO NOT** twist hand grips as they are not meant to be twisted. Failure to do this could result in a fall.
- **DO NOT** hang anything from the frame of the rolling walker other than the basket included or the tote bag accessory. The basket has a weight limitation of 11 lbs (5 kg). The tote bag has a weight limitation of 10 lbs (4.5 kg). Items should be placed completely in the basket or tote bag and should not protrude. This may cause the rollator to tip, resulting in bodily injury or rollator damage.
- The backrest is intended to provide back support as you sit on the seat. The backrest is not intended to support full weight. Avoid leaning while sitting as this may result in fall.

Assembly Instructions

IMPORTANT: Read all operating and safety instructions before using your rollator.

Note: You may wish to request a set up and operating demonstration with your local Medical Equipment provider or pharmacy.

Figure 1: Assembled Rollator


Figure 2:
Knob Screw

Rollator Components

- 1 Rollator frame
- 2 Front wheels
- 2 Handlebars
- 2 Rear wheels
- 2 Knob screws
- 1 Backrest with push-buttons
- 2 Round hand screws
- 1 Basket with handle

Assembled Rollator

- A. Handgrips
- B. Adjustable Handlebar
- C. Padded Backrest
- D. Brake Cables
- E. Padded Seat
- F. Support Rod
- G. 8" (203 mm) Front Wheels
- H. 8" (203 mm) Rear Wheels
- I. Folding Bracket
- J. C Clamps
- K. Knob Screw

Assembly Instructions

1. Carefully remove the rollator components from the shipping carton. Refer to diagram **(Figure 1)** as needed.
2. Examine all parts for shipping damage. **DO NOT** use the rollator if there is apparent damage or if any parts (other than the basket) are missing. If there are missing or damaged parts, please call Medline Customer Care toll free at 1-800-MEDLINE.
3. Lay the rollator down so that the seat is facing the floor. Cut the plastic ties from the rear wheels and unfold the legs. Straighten the legs until the rollator unfolds.
4. Flip the rollator over so that the seat faces up. Insert front wheel (G) into frame. Insert C Clamp (J) into frame.
5. Unfold the rollator by pushing down on the support rod (F) to spread apart the frame until the folding brackets (I) are in full extension. Push the seat down into place.
6. Insert the handlebars (B) so that the handgrips (A) are aligned with the rear wheels. The brake cables should be on the outside of the handlebars and frame. Ensure excess brake cable loops out the front of the rollator.
7. Adjust the height of the handlebars (B) for proper fit. To determine proper height, stand upright behind the rollator with shoulders relaxed, arms hanging at your side and eyes looking straight ahead. Standing in this position, the handgrips (A) of the rollator should be adjusted to a height even with the wrists. This will position the arms at an approximate 20° to 30° bend when using the rollator. Secure the handlebars at the closest hold to the desired height. Insert and tighten the triangular hand screws until secure. **See Figure 2. (Important: Hand Tighten Only!)**
8. Insert the padded backrest (C) so that the padded portion faces the front of the unit. Push down on the buttons on the side of the backrest until the spring pins lock into place.
9. Place the basket under the seat.

Operating Instructions

Opening the Rollator

Unfold the Medline rollator by pushing down on the support rod **(F)** to spread the frame apart until the folding brackets **(I)** are in full extension. Push the seat down into place.

Folding the Rollator

Remove the basket from under the seat. Pull up on the folding bracket tie-bar on the support rod **(F)** until the frame is fully collapsed.

Checking the Brakes

- Before each use, check to make sure the brakes are working properly.
- Applying pressure on brake handles using a pulling up motion the brakes will stop the wheels momentarily.
- Pushing down on the brake handles will lock the wheels in place permanently.

Adjusting the Brakes

1. Loosen the cable nut.
2. While holding the brake to the wheel.
3. Pull down on the brake cable to lengthen the brake cord.
4. Retighten the nut. Retest and repeat as necessary for good brake function.

Caution: After making any adjustments to the rollator brakes, you must test the brakes before the rollator is used. See Maintenance instructions.

Operating Instructions

Using the rollator as a walker:

For proper positioning and use of your rollator, imagine a line drawn across the back of the rear wheels. Part of one foot should be slightly ahead of that line as you walk.

Using the rollator as a seat:

- Push the rollator into the position where you want to sit
- **LOCK** the brakes **BEFORE** sitting.

To Operate the Brakes:

- To slow down, squeeze the brake handles up.
- To lock the brakes, push down on the brake handles until you feel or hear the brakes click into locked position. If the brakes have been adjusted but the rear wheels on the rollator move while the brakes are in the locked position, **DO NOT** use the rollator. Please call our toll-free customer service number at 1-800-MEDLINE for returns or repairs.
- To release the brakes from the locked position, simply pull up on the brake handles.

Safety Warnings: Also See Instructions

Before using the seat portion of your rollator, **ALWAYS** lock the brakes **FIRST**. **DO NOT** navigate the rollator while sitting on the seat. The brakes should always be locked when the seat is being used. Risk of fall and serious injury may result if brakes are **NOT** locked while the user is sitting on the seat.

Maintenance

Your rollator should be checked periodically to ensure the brakes are functioning properly and that all nuts and bolts are secure.

Periodically test the brakes by both walking and squeezing the brakes and by stopping, pushing down on the brakes to a locked position, and then trying to push the walker.

If the brakes do not lock the wheels, the brakes need to be adjusted. If adjustment does not tighten the brakes, **DO NOT** use the rollator.

Keep your receipt!

Contact us at 1-800-MEDLINE, Monday – Friday 8am to 5pm CST if you encounter any problems or for information regarding repairs or maintenance.

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LIMITED WARRANTY

Your Medline product is warranted to be free of defects in materials and workmanship for the duration of use by the original purchaser. This device was built to precise standards and carefully inspected prior to shipment. This limited warranty is an expression of our confidence in the materials and workmanship of our products and assures years of convenient, dependable service. In the event of a defect covered by this warranty, we will, at our option, repair or replace the device. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters and grips, which are subject to normal wear and replacement.

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Medline Canada
1-800-396-6996
www.medline.ca | Canada@medline.com

Medline Mexico
01-800-831-0898
www.medlinemexico.com | mexico@medline.com

Some products may not be available for sale in Mexico or Canada. Please contact your Account Representative for additional information. We reserve the right to correct any errors that may occur within this manual. ©2010 Medline industries, Inc. Guardian is a registered trademark of Medline Industries, Inc.

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