

A system you can rely on

Philips Respironics InnoSpire Essence nebulizer compressor system

A reliable choice

When reliability is key, choose InnoSpire Essence. Essence, part of our InnoSpire family of compressor nebulizer systems, provides fast, efficient and reliable drug delivery to patients at home.

Proven nebulizer technology

Essence is matched with our SideStream disposable nebulizer to deliver a high quality aerosol in a short treatment time. With more than 20 years of aerosol therapy excellence, you can rest assured that patients will be provided with a high quality aerosol every time.

Key features

- Partnered with proven
 SideStream technology
- Short treatment times in 6 8 minutes*
- Compact design





Reliable aerosol delivery you can trust

Essence is combined with our highly efficient SideStream nebulizer. SideStream has been designed to draw in more air, via an active venturi system, resulting in faster drug delivery and shorter treatment times. It also has a unique diamond jet design that helps to minimize wear from frequent use and aids consistent delivery of a high quality aerosol.



Performance data	
MMAD	2.9 µm
Fine particle fraction	77% below 5 μm
Treatment time (2.5 ml charge)	6 - 8 mins

Product description (with voltage/frequency)	
Essence 230 V/50 Hz EU plug	1099967
Essence 220 V/60 Hz Korean plug	1099965
Essence 220 V/60 Hz UK plug	1099964
Essence 230 V/50 Hz UK plug	1112261
Essence 220 V/60 Hz US plug	1102174
SideStream reusable year pack	2010A
SideStream reusable kit with adult mask	1223A
SideStream reusable kit with pediatric mask	1224A
SideStream reusable kit with mouthpiece	1225A
Air inlet filters (pack of 4)	1102088
SideStream disposable one month pack (pack of 40)	2350
SideStream disposable with adult mask (pack of 50)	4446
SideStream disposable with pediatric mask (pack of 50)	4447
SideStream disposable with mouthpiece (pack of 50)	4448

230 V/50 Hz, 1 Amp	220 V/60 Hz, 1.2 Amp
317 kPa (46 psi)	296 kPa (43 psi)
6 lpm @ 10 psig	7 lpm @ 10 psig
8.0 lpm	9.3 lpm
Thermal fuse; functioning temperature 150°C	
1.5 kg	
165 x 165 x 108 mm	
58 ±3 dBA	
	317 kPa (46 psi) 6 lpm @ 10 psig 8.0 lpm Thermal fuse; function 1.5 kg 165 x 165 x 108 mm

^{*}data on file at Philips Respironics



Respironics, InnoSpire, Essence and SideStream are trademarks of Koninklijke Philips N.V. and its affiliates, all rights reserved.

Scan with your smartphone or tablet to go directly to our website.

Please visit www.philips.com/innospire for more information



© 2014 Koninklijke Philips N.V.All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

BROUDY HT 5/17/14 MCI 4106223 PN 1103686

Philips Healthcare is part of Royal Philips

Philips Respironics Asia Pacific +65 6882 5282 Philips Respironics Australia +61 (2) 9947 0440 Philips Respironics United Kingdom +44 800 1300 845 Respironics Respiratory Drug Delivery (UK) Ltd Chichester Business Park, City Fields Way Tangmere, Chichester P020 2FT UK

Customer Service +1 724 387 4000 www.philips.com/respironics